


DISEASE EXOME, A POWERFUL DIAGNOSTIC TOOL: POST-MORTEM DIAGNOSIS OF DYSKERATOSIS CONGENITA



Pinto-Basto J, Melo F, Almeida LS, Garrido A, Barbosa C, Pereira SM, Cerqueira R, Tavares P



Lígia S Almeida

ligia.almeida@cgcgenetics.com

Senior Clinical Molecular Geneticist,
Molecular Diagnostics & Clinical Genomics Laboratories, CGC Genetics

19ª Reunião Anual SPGH

7 de Novembro 2015

DYSKERATOSIS CONGENITA, DKC

> BACKGROUND

1. Rare ectodermal dysplasia that often presents with the classic triad of:

- Nail dysplasia,
- Skin pigmentation changes
- Oral leukoplakia associated with a high risk of bone marrow failure (BMF) and cancer.

2. Onset and progression of manifestations of DKC varies:



3. Diagnosis:

- Telomere length - cytogenetics analysis (ex: FISH)
- **Molecular Diagnosis** - mutations known to cause DKC in:

CTC1

DKC1

TERC

TERT

TINF2

NHP2

NOP10

WRAP53

There are several variants of the disorder, each associated with a different pattern of inheritance

CASE REPORT

> OVERVIEW

PATIENT:


- > Deceased 60 year-old male
- > Undiagnosed systemic disease (age at onset - 50 years):
 - liver cirrhosis
 - pulmonary fibrosis
 - sick sinus syndrome
 - thrombocytopenia
- > Multiple tests had been previously performed with **no** diagnosis.

DNA from *post-mortem* tissue was obtained and *DISEASE EXOME* performed.

DISEASE EXOME

> DISEASE EXOME BY CGC GENETICS

- Uncertain or complex phenotype
- Unclear diagnosis




DISEASE EXOME panel is based on:

1. Human Gene Mutation Database (HGMD) (www.hgmd.cf.ac.uk/ac/index.php)
2. Online Mendelian Inheritance in Man (OMIM) (www.omim.org)
3. GeneTests (www.genetests.org)
4. Illumina TruSight Sequencing Panels (www.illumina.com/trusight)
5. Other commercially available sequencing panels


METHODOLOGY PROCESS

> NEXT GENERATION SEQUENCING (NGS)



METHODOLOGY PROCESS

> NEXT GENERATION SEQUENCING (NGS)




> Capture of target regions using oligonucleotide probes
TruSightOne, Illumina

> NGS of a panel composed by 4,813 clinically-relevant genes
(Disease Exome)
MiSeq, Illumina



METHODOLOGY PROCESS

> NEXT GENERATION SEQUENCING (NGS)




- > Alignment and variant calling
BWA and GATK
- > Variants with MAF<1%
Bioinformatic analysis tools

Alignment example:

Reads	ATGGCATTGCAATTTGACAT
	TGGCATTGCAATTTG
	AGATGGTATTG
	GATGGCATTGCAA
	GCATTGCAATTTGAC
	ATGGCATTGCAATT
AGATGGCATTGCAATTTG	
Reference Genome	AGATGGTATTGCAATTTGACAT



METHODOLOGY PROCESS

> NEXT GENERATION SEQUENCING (NGS)




METHODOLOGY PROCESS

> NEXT GENERATION SEQUENCING (NGS)



METHODOLOGY PROCESS

> NEXT GENERATION SEQUENCING (NGS)



GENETIC ANALYSIS


> PRIORITIZATION OF DISEASE EXOME DATA

Exclusion of:

- low quality reads
- non-genic, intronic and synonymous variants
- variants in in-house variant database

Inclusion of:

- Variants with less than 1%
- Variants in exonic, splice-site and canonical splice site regions



GENETIC ANALYSIS

> DISEASE EXOME BY CGC GENETICS

Heterozygous variant in *TERT* gene
NM_198253.2: c.1492G>A (p.Gly498Arg)

- > Variant not yet described in the literature
- Not detected:** dbSNP, ExAC, 1KGenomes, ESP
- > Affects highly conserved residue
- > Bioinformatic analysis - likely pathogenic variant
 - **PoplyPhen-2:** probably damaging
 - **SIFT:** deleterious

TERT {
> **Autosomal Dominant Dyskeratosis Congenita**
(type 2; MIM 613989)
> **Autosomal Recessive Dyskeratosis Congenita**
(type 4; MIM 613989)


GENETIC ANALYSIS

> TERT - Telomerase reverse transcriptase (Chr 5)




NM_198253.2:
c.1492G>A (p.Gly498Arg)

> NGS



> SANGER CONFIRMATION



CONCLUSION

- ***post-mortem* diagnosis of autosomal dominant DKC type-2**
- **Novel variant in *TERT* gene, expanding the mutational spectrum of *TERT*-related DKC**
- **Highlights the importance of molecular diagnosis, even *post-mortem*, as establishing the molecular etiology allows proper genetic counselling to at-risk relatives**
- **Clinical exome sequencing applied in solving complex diagnostic cases**
 - **Directed to clinically relevant genes**
 - **High diagnostic yield**
 - **Clinical integration and interpretation**
 - **Lower cost and more useful comparatively to other exome testing**



ACKNOWLEDGMENTS

Purificação Tavares (CEO)
Molecular Diagnostics & Clinical Genomics Laboratories

Jorge Pinto Basto

Rita Cerqueira

Marisa Teixeira

Fátima Torres

Isabel Miguel

Ana Oliveira

Alexandra Lopes

Lígia Lameiras

Catarina Ribeiro

Luena Pitrez

Sílvia Pereira

Inês Castro

Liliana Castro

Sofia Lisboa

Sofia Pereira

Daniela Teixeira

NGS Team

Filipa Melo

Luísa Bastos

Andreia Pinto

Mariana Ferreira

Andreia Pinto

Viviana Silva

Susana Neto

Marisa Paulino

Centro Hospitalar Tondela-Viseu, EPE

António José Garrido

Cristina Barbosa

Thank you for your
Attention!

